




VCP Antibody

Product Code	CSB-PA060060
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55072
Immunogen	Synthesized peptide derived from Human VCP around the non-phosphorylation site of S352.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	VCP; Transitional endoplasmic reticulum ATPase; TER ATPase; 15S Mg(2+)-ATPase p97 subunit; Valosin-containing protein; VCP
Immunogen Species	Homo sapiens (Human)
Target Names	VCP

Image



Western Blot analysis of K562 cells using VCP Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.