



Phospho-VHL (S68) Antibody

Product Code	CSB-PA060066
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40337
Immunogen	Synthesized peptide derived from Human VHL around the phosphorylation site of S68.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	VHL; Von Hippel-Lindau disease tumor suppressor; Protein G7; pVHL
Immunogen Species	Homo sapiens (Human)
Target Names	VHL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.