





EZR Antibody

Product Code	CSB-PA060072
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15311
Immunogen	Synthesized peptide derived from Human Ezrin around the non-phosphorylation site of Y146.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	EZR; VIL2; Ezrin; Cytovillin; Villin-2; p81
Immunogen Species	Homo sapiens (Human)
Target Names	EZR

Image



Western Blot analysis of Jurkat cells using Ezrin Polyclonal Antibody



Western Blot analysis of Jurkat cells using Ezrin Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.