



Phospho-WAS (Y290) Antibody

Product Code	CSB-PA060081
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42768
Immunogen	Synthesized peptide derived from Human WASP around the phosphorylation site of Y290.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	WAS; IMD2; Wiskott-Aldrich syndrome protein; WASp
Immunogen Species	Homo sapiens (Human)
Target Names	WAS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.