



YWHAZ/YWHAB Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA060095 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P63104/P31946 |
| Immunogen | Synthesized peptide derived from Human 14-3-3 β/ζ around the non-phosphorylation site of S184/186. |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse, Rat |
| Tested Applications | WB, ELISA; WB:1:500-1:2000, ELISA:1:20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | YWHAZ; 14-3-3 protein zeta/delta; Protein kinase C inhibitor protein 1; KCIP-1; YWHAB; 14-3-3 protein beta/alpha; Protein 1054; Protein kinase C inhibitor protein 1; KCIP-1 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | YWHAZ/YWHAB |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |