




# Phospho-SF1 (S82) Antibody

<b>Product Code</b>	CSB-PA060104
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q15637
<b>Immunogen</b>	Synthesized peptide derived from Human Splicing factor 1 around the phosphorylation site of S82.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Monkey
<b>Tested Applications</b>	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	SF1; ZFM1; ZNF162; Splicing factor 1; Mammalian branch point-binding protein; BBP; mBBP; Transcription factor ZFM1; Zinc finger gene in MEN1 locus; Zinc finger protein 162
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	SF1

## Image



Western blot analysis of SH-SY5Y 293T using p-Splicing factor 1 (S82) antibody.