



ZNF7 Antibody

Product Code	CSB-PA060108
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17097
Immunogen	Synthesized peptide derived from the N-terminal region of Human KOX4.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ZNF7; KOX4; Zinc finger protein 7; Zinc finger protein HF.16; Zinc finger protein KOX4
Immunogen Species	Homo sapiens (Human)
Target Names	ZNF7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.