




CA2 Antibody

Product Code	CSB-PA060111
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00918
Immunogen	Synthesized peptide derived from the C-terminal region of Human CA II.
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CA2; Carbonic anhydrase 2; Carbonate dehydratase II; Carbonic anhydrase C; CAC; Carbonic anhydrase II; CA-II
Immunogen Species	Homo sapiens (Human)
Target Names	CA2

Image



Western blot analysis of Mouse-kidney mouse-brain 293T hela lysis using CA II antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.