



ZNF76 Antibody

Product Code	CSB-PA060113
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36508
Immunogen	Synthesized peptide derived from the Internal region of Human ZNF76.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ZNF76; D6S229E; ZNF523; Zinc finger protein 76; Zinc finger protein 523
Immunogen Species	Homo sapiens (Human)
Target Names	ZNF76
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.