



Phospho-TUBA1A/TUBA1B/TUBA1C/TUBA3C/TUBA3E/TUBA4A (Y272) Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA060125 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q71U36/P68363/Q9BQE3/Q13748/Q6PEY2/P68366 |
| Immunogen | Synthesized peptide derived from Human Tubulin α 1/3/4 around the phosphorylation site of Y272. |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse, Rat |
| Tested Applications | WB, ELISA; WB:1:500-1:2000, ELISA:1:40000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | TUBA1A; TUBA3; Tubulin alpha-1A chain; Alpha-tubulin 3; Tubulin B-alpha-1; Tubulin alpha-3 chain; TUBA1B; Tubulin alpha-1B chain; Alpha-tubulin ubiquitous; Tubulin K-alpha-1; Tubulin alpha-ubiquitous chain; TUBA1C; TUBA6; Tubulin alpha-1C c |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | TUBA1A/TUBA1B/TUBA1C/TUBA3C/TUBA3E/TUBA4A |