




SLC52A2 Antibody

Product Code	CSB-PA060150
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HAB3
Immunogen	Synthesized peptide derived from the Internal region of Human GPR172A.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SLC52A2; GPR172A; PAR1; RFT3; Solute carrier family 52; riboflavin transporter, member 2; Porcine endogenous retrovirus A receptor 1; PERV-A receptor 1; Protein GPR172A; Riboflavin transporter 3; hRFT3
Immunogen Species	Homo sapiens (Human)
Target Names	SLC52A2

Image



Western Blot analysis of A549 cells using GPR172A Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.