



CNTNAP3 Antibody

Product Code	CSB-PA060170
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZ76
Immunogen	Synthesized peptide derived from the C-terminal region of Human CASPR3.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CNTNAP3; CASPR3; KIAA1714; Contactin-associated protein-like 3; Cell recognition molecule Caspr3
Immunogen Species	Homo sapiens (Human)
Target Names	CNTNAP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.