



CSRNP3 Antibody

Product Code	CSB-PA060177
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WYN3
Immunogen	Synthesized peptide derived from the Internal region of Human CSRNP3.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Monkey
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CSRNP3; FAM130A2; TAIP2; Cysteine/serine-rich nuclear protein 3; CSRNP-3; Protein FAM130A2; TGF-beta-induced apoptosis protein 2; TAIP-2
Immunogen Species	Homo sapiens (Human)
Target Names	CSRNP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.