



NDFIP1 Antibody

Product Code	CSB-PA060189
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BT67
Immunogen	Synthesized peptide derived from the C-terminal region of Human Ndfip1.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NDFIP1; N4WBP5; PSEC0192; PSEC0223; NEDD4 family-interacting protein 1; Breast cancer-associated protein SGA-1M; NEDD4 WW domain-binding protein 5; Putative MAPK-activating protein PM13; Putative NF-kappa-B-activating protein 164; Putative
Immunogen Species	Homo sapiens (Human)
Target Names	NDFIP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.