



OR11H1/OR11H2/OR11H12 Antibody

Product Code	CSB-PA060197
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NG94/Q8NH07/B2RN74/Q6IEX0
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 11H1/11H2/11H12.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR11H1; Olfactory receptor 11H1; Olfactory receptor OR22-1; OR11H2; C14orf15; OR11H2P; Olfactory receptor 11H2; Olfactory receptor OR14-1; OR11H12; Olfactory receptor 11H12Q6IEX0
Immunogen Species	Homo sapiens (Human)
Target Names	OR11H1/OR11H2/OR11H12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.