



OR4F17/OR4F4/OR4F5 Antibody

Product Code	CSB-PA060198
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGA8/Q96R69/Q8NH21
Immunogen	Synthesized peptide derived from the Internal region of Human Olfactory receptor 4F4/4F5/4F17.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR4F17; OR4F11P; OR4F18; OR4F19; Olfactory receptor 4F17; Olfactory receptor 4F11; Olfactory receptor 4F18; Olfactory receptor 4F19; OR4F4; Olfactory receptor 4F4; HS14a-1-A; Olfactory receptor OR19-3; OR4F5; Olfactory receptor 4F5
Immunogen Species	Homo sapiens (Human)
Target Names	OR4F17/OR4F4/OR4F5