



# OR4F17/OR4F4/OR4F5 Antibody

<b>Product Code</b>	CSB-PA060198
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8NGA8/Q96R69/Q8NH21
<b>Immunogen</b>	Synthesized peptide derived from the Internal region of Human Olfactory receptor 4F4/4F5/4F17.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	OR4F17; OR4F11P; OR4F18; OR4F19; Olfactory receptor 4F17; Olfactory receptor 4F11; Olfactory receptor 4F18; Olfactory receptor 4F19; OR4F4; Olfactory receptor 4F4; HS14a-1-A; Olfactory receptor OR19-3; OR4F5; Olfactory receptor 4F5
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	OR4F17/OR4F4/OR4F5
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.