



FZD4 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA060236 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9ULV1 |
| Immunogen | Synthesized peptide derived from the Internal region of Human Frizzled-4. |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse, Rat |
| Tested Applications | IF, ELISA; IF:1:200-1:1000, ELISA:1:10000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | FZD4; Frizzled-4; Fz-4; hFz4; FzE4; CD antigen CD344 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | FZD4 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |