



Histone H2A Antibody

Product Code	CSB-PA060244
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C0S8/P04908/P20671/Q96KK5/Q6FI13/Q7L7L0
Immunogen	Synthesized peptide derived from Human Histone H2A around the non-phosphorylation site of T121.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HIST1H2AG; H2AFP; HIST1H2AI; H2AFC; HIST1H2AK; H2AFD; HIST1H2AL; H2AFI; HIST1H2AM; H2AFN; Histone H2A type 1; H2A.1; Histone H2A/p; HIST1H2AB; H2AFM; HIST1H2AE; H2AFA; Histone H2A type 1-B/E; Histone H2A.2; Histone H2A/a; Histone H2A/m; HIS
Immunogen Species	Homo sapiens (Human)
Target Names	HIST1H2AG/HIST1H2AB/HIST1H2AD/HIST1H2AH/HIST2H2AA3/HIST3H2A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.