



ELOVL3 Antibody

Product Code	CSB-PA060250
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HB03
Immunogen	Synthesized peptide derived from the Internal region of Human ELOVL3.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ELOVL3; CIG30; Elongation of very long chain fatty acids protein 3; 3-keto acyl-CoA synthase ELOVL3; Cold-inducible glycoprotein of 30 kDa; ELOVL fatty acid elongase 3; ELOVL FA elongase 3
Immunogen Species	Homo sapiens (Human)
Target Names	ELOVL3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.