



# Histone H3 Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA060259   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Uniprot No.</b>         | P68431/Q71DI3/P84243   |
| <b>Immunogen</b>           | Synthesized peptide derived from Human Histone H3 around the non-phosphorylation site of T11.  |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human, Mouse, Rat  |
| <b>Tested Applications</b> | IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000  |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Non-conjugated   |
| <b>Storage Buffer</b>      | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.  |
| <b>Purification Method</b> | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.  |
| <b>Isotype</b>             | IgG  |
| <b>Alias</b>               | HIST1H3A; H3FA; HIST1H3B; H3FL; HIST1H3C; H3FC; HIST1H3D; H3FB; HIST1H3E; H3FD; HIST1H3F; H3FI; HIST1H3G; H3FH; HIST1H3H; H3FK; HIST1H3I; H3FF; HIST1H3J; H3FJ; Histone H3.1; Histone H3/a; Histone H3/b; Histone H3/c; Histone H3/d; Histone H3 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Target Names</b>        | HIST1H3A/HIST2H3A/H3F3A  |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.  |