



# Cleaved-CASP4/CASP5 (D270/D311) Antibody

<b>Product Code</b>	CSB-PA060272
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P49662/P51878
<b>Immunogen</b>	Synthesized peptide derived from the Internal region of Human Caspase-4/5 p20.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	CASP4; ICH2; Caspase-4; CASP-4; ICE(rel)-II; Protease ICH-2; Protease TX; CASP5; ICH3; Caspase-5; CASP-5; ICE(rel)-III; Protease ICH-3; Protease TY
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	CASP4/CASP5
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.