



MAD1L1 Antibody

Product Code	CSB-PA060283
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6D9
Immunogen	Synthesized peptide derived from Human MAD1 around the non-phosphorylation site of S428.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MAD1L1; MAD1; TXBP181; Mitotic spindle assembly checkpoint protein MAD1; Mitotic arrest deficient 1-like protein 1; MAD1-like protein 1; Mitotic checkpoint MAD1 protein homolog; HsMAD1; hMAD1; Tax-binding protein 181
Immunogen Species	Homo sapiens (Human)
Target Names	MAD1L1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.