



FGFBP2 Antibody

Product Code	CSB-PA060288
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYJ0
Immunogen	Synthesized peptide derived from the C-terminal region of Human FGF-BP2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FGFBP2; KSP37; Fibroblast growth factor-binding protein 2; FGF-BP2; FGF-binding protein 2; FGFBP-2; 37 kDa killer-specific secretory protein; Ksp37; HBp17-related protein; HBp17-RP
Immunogen Species	Homo sapiens (Human)
Target Names	FGFBP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.