




Phospho-CASP6 (S257) Antibody

Product Code	CSB-PA060289
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55212
Immunogen	Synthesized peptide derived from Human Caspase-6 around the phosphorylation site of S257.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CASP6; MCH2; Caspase-6; CASP-6; Apoptotic protease Mch-2
Immunogen Species	Homo sapiens (Human)
Target Names	CASP6

Image



Western Blot analysis of Jurkat cells using Phospho-Caspase-6 (S257) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.