



PLA2G6 Antibody

Product Code	CSB-PA060291
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60733
Immunogen	Synthesized peptide derived from the Internal region of Human Group VI iPLA2.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PLA2G6; PLPLA9; 85/88 kDa calcium-independent phospholipase A2; Cal-PLA2; Group VI phospholipase A2; GVI PLA2; Intracellular membrane-associated calcium-independent phospholipase A2 beta; iPLA2-beta; Patatin-like phospholipase domain-contai
Immunogen Species	Homo sapiens (Human)
Target Names	PLA2G6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.