



NCK2 Antibody

Product Code	CSB-PA070023
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43639
Immunogen	Synthesized peptide derived from the C-terminal region of Human Nck-2.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NCK2; GRB4; Cytoplasmic protein NCK2; Growth factor receptor-bound protein 4; NCK adaptor protein 2; Nck-2; SH2/SH3 adaptor protein NCK-beta
Immunogen Species	Homo sapiens (Human)
Target Names	NCK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.