









COX4I2 Antibody

Product Code	CSB-PA070044
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96KJ9
Immunogen	Synthesized peptide derived from the Internal region of Human COX4I2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Conjugate Storage Buffer	Non-conjugated Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. The antibody was affinity-purified from rabbit antiserum by affinity-
Storage Buffer Purification Method	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Storage Buffer Purification Method Isotype	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. IgG COX4I2; COX4L2; Cytochrome c oxidase subunit 4 isoform 2; mitochondrial;