



ABHD14B Antibody

Product Code	CSB-PA070049
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96IU4
Immunogen	Synthesized peptide derived from the C-terminal region of Human ABHD14B.
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ABHD14B; CIB; Abhydrolase domain-containing protein 14B; CCG1-interacting factor B
Immunogen Species	Homo sapiens (Human)
Target Names	ABHD14B