



# Phospho-IKBKG (S85) Antibody

<b>Product Code</b>	CSB-PA070068
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9Y6K9
<b>Immunogen</b>	Synthesized peptide derived from Human IKK $\gamma$ around the phosphorylation site of S85.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	IKBKG; FIP3; NEMO; NF-kappa-B essential modulator; NEMO; FIP-3; I $\kappa$ B kinase-associated protein 1; IKKAP1; Inhibitor of nuclear factor kappa-B kinase subunit gamma; I-kappa-B kinase subunit gamma; IKK-gamma; IKKG; I $\kappa$ B kinase subunit gamma; NF
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	IKBKG
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.