




Phospho-IKBKG (S31) Antibody

Product Code	CSB-PA070069
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6K9
Immunogen	Synthesized peptide derived from Human IKK γ around the phosphorylation site of S31.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	IKBKG; FIP3; NEMO; NF-kappa-B essential modulator; NEMO; FIP-3; I κ B kinase-associated protein 1; IKKAP1; Inhibitor of nuclear factor kappa-B kinase subunit gamma; I-kappa-B kinase subunit gamma; IKK-gamma; IKKG; I κ B kinase subunit gamma; NF
Immunogen Species	Homo sapiens (Human)
Target Names	IKBKG

Image



Western Blot analysis of 293 cells using Phospho-IKK γ (S31) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.