



CDC14A Antibody

Product Code	CSB-PA070087
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UNH5
Immunogen	Synthesized peptide derived from the Internal region of Human Cdc14a phosphatase.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CDC14A; Dual specificity protein phosphatase CDC14A; CDC14 cell division cycle 14 homolog A
Immunogen Species	Homo sapiens (Human)
Target Names	CDC14A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.