



PLA2G4C Antibody

Product Code	CSB-PA070099
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UP65
Immunogen	Synthesized peptide derived from the Internal region of Human cPLA2-γ.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PLA2G4C; Cytosolic phospholipase A2 gamma; cPLA2-gamma; Phospholipase A2 group IVC
Immunogen Species	Homo sapiens (Human)
Target Names	PLA2G4C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.