



CBFA2T3 Antibody

Product Code	CSB-PA070109
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75081
Immunogen	Synthesized peptide derived from the Internal region of Human ETO-2.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CBFA2T3; MTG16; MTGR2; ZMYND4; Protein CBFA2T3; MTG8-related protein 2; Myeloid translocation gene on chromosome 16 protein; hMTG16; Zinc finger MYND domain-containing protein 4
Immunogen Species	Homo sapiens (Human)
Target Names	CBFA2T3