



# Recombinant Mouse Mannose-binding protein C (Mbl2)

<b>Product Code</b>	CSB-YP013540MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P41317
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	ET LTEGVQNSCP VVTCSSPGLN GFPGKDGRDG AKGEKGEPGQ GLRGLQGPPG KVGPTGPPGN PGLKGAVGPK GDRGDRAEFD TSEIDSEIAA LRSELRALRN WVLFSLSEKV GK KYFVSSVK KMSLDRVKAL CSEFQGSVAT PRNAEENSAI QKVAKDIAYL GITDVRVEGS FEDLTGNRVR YTNWNDGEPN NTGDGEDCVV ILGNGKWNDV PCSDSFLAIC EFS D
<b>Source</b>	Yeast
<b>Target Names</b>	Mbl2
<b>Protein Names</b>	Recommended name: Mannose-binding protein C Short name= MBP-C Alternative name(s): Mannan-binding protein RA-reactive factor P28A subunit Short name= RARF/P28A
<b>Expression Region</b>	19-244
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes the soluble mannose-binding lectin or mannose-binding protein found in serum. The protein encoded belongs to the collectin family and is an important element in the innate immune system. The protein recognizes mannose and N-acetylglucosamine on many microorganisms, and is capable of activating the classical complement pathway. Deficiencies of this gene have been associated with susceptibility to autoimmune and infectious diseases.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.