



# B3GALT2 Antibody

<b>Product Code</b>	CSB-PA070132
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O43825
<b>Immunogen</b>	Synthesized peptide derived from the C-terminal region of Human $\beta$ -1,3-Gal-T2.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Tested Applications</b>	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:40000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	B3GALT2; Beta-1; 3-galactosyltransferase 2; Beta-1,3-GalTase 2; Beta3Gal-T2; Beta3GalT2; UDP-galactose:2-acetamido-2-deoxy-D-glucose 3beta-galactosyltransferase 2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	B3GALT2
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.