




Phospho-FADD (S194) Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA070144 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q13158 |
| Immunogen | Synthesized peptide derived from Human FADD around the phosphorylation site of S194. |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | FADD; MORT1; GIG3; Protein FADD; FAS-associated death domain protein; FAS-associating death domain-containing protein; Growth-inhibiting gene 3 protein; Mediator of receptor induced toxicity |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | FADD |

Image



Western blot analysis of MOUSE-KIDNEY
MOUSE-BRAIN MOUSE-SPLEEN MOUSE-
LUNG using p-FADD (S194) antibody.