




Phospho-FADD (S194) Antibody

Product Code	CSB-PA070144
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13158
Immunogen	Synthesized peptide derived from Human FADD around the phosphorylation site of S194.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FADD; MORT1; GIG3; Protein FADD; FAS-associated death domain protein; FAS-associating death domain-containing protein; Growth-inhibiting gene 3 protein; Mediator of receptor induced toxicity
Immunogen Species	Homo sapiens (Human)
Target Names	FADD

Image



Western blot analysis of MOUSE-KIDNEY, MOUSE-BRAIN, MOUSE-SPLEEN, and MOUSE-LUNG using p-FADD (S194) antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.