



ALDH1A2 Antibody

Product Code	CSB-PA070161
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94788
Immunogen	Synthesized peptide derived from the C-terminal region of Human ALDH1A2.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IHC:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ALDH1A2; RALDH2; Retinal dehydrogenase 2; RALDH 2; RaIDH2; Aldehyde dehydrogenase family 1 member A2; Retinaldehyde-specific dehydrogenase type 2; RALDH(II)
Immunogen Species	Homo sapiens (Human)
Target Names	ALDH1A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.