




Phospho-WASF1 (Y125) Antibody

Product Code	CSB-PA070180
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92558
Immunogen	Synthesized peptide derived from Human WAVE1 around the phosphorylation site of Y125.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	WASF1; KIAA0269; SCAR1; WAVE1; Wiskott-Aldrich syndrome protein family member 1; WASP family protein member 1; Protein WAVE-1; Verprolin homology domain-containing protein 1
Immunogen Species	Homo sapiens (Human)
Target Names	WASF1

Image



Western Blot analysis of 3T3 cells using Phospho-WAVE1 (Y125) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.