



CCND2 Antibody

Product Code	CSB-PA070182
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30279
Immunogen	Synthesized peptide derived from Human Cyclin D2 around the non-phosphorylation site of T280.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CCND2; G1/S-specific cyclin-D2
Immunogen Species	Homo sapiens (Human)
Target Names	CCND2