



Phospho-CCNE1 (T77) Antibody

Product Code	CSB-PA070188
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24864
Immunogen	Synthesized peptide derived from Human Cyclin E1 around the phosphorylation site of T77.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CCNE1; CCNE; G1/S-specific cyclin-E1
Immunogen Species	Homo sapiens (Human)
Target Names	CCNE1