




Phospho-STAM2 (Y192) Antibody

Product Code	CSB-PA070192
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75886
Immunogen	Synthesized peptide derived from Human STAM2 around the phosphorylation site of Y192.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	STAM2; HBP; Signal transducing adapter molecule 2; STAM-2; Hrs-binding protein
Immunogen Species	Homo sapiens (Human)
Target Names	STAM2

Image



Western Blot analysis of 3T3 cells using Phospho-STAM2 (Y192) Polyclonal Antibody