




Phospho-STAM2 (Y192) Antibody

Product Code	CSB-PA070192
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75886
Immunogen	Synthesized peptide derived from Human STAM2 around the phosphorylation site of Y192.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	STAM2; HBP; Signal transducing adapter molecule 2; STAM-2; Hrs-binding protein
Immunogen Species	Homo sapiens (Human)
Target Names	STAM2

Image



Western Blot analysis of 3T3 cells using Phospho-STAM2 (Y192) Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.