




Phospho-RPS6KA4 (T568) Antibody

Product Code	CSB-PA070193
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75676
Immunogen	Synthesized peptide derived from Human MSK2 around the phosphorylation site of T568.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RPS6KA4; MSK2; Ribosomal protein S6 kinase alpha-4; S6K-alpha-4; 90 kDa ribosomal protein S6 kinase 4; Nuclear mitogen- and stress-activated protein kinase 2; Ribosomal protein kinase B; RSKB
Immunogen Species	Homo sapiens (Human)
Target Names	RPS6KA4

Image



Western Blot analysis of 293 cells using Phospho-MSK2 (T568) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.