



Phospho-ZNF598 (Y306) Antibody

Product Code	CSB-PA070224
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UK7
Immunogen	Synthesized peptide derived from Human ZFP598 around the phosphorylation site of Y306.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ZNF598; Zinc finger protein 598
Immunogen Species	Homo sapiens (Human)
Target Names	ZNF598
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.