



Phospho-CCNE2 (T392) Antibody

Product Code	CSB-PA070234
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O96020
Immunogen	Synthesized peptide derived from Human Cyclin E2 around the phosphorylation site of T392.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CCNE2; G1/S-specific cyclin-E2
Immunogen Species	Homo sapiens (Human)
Target Names	CCNE2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.