



NEK9 Antibody

Product Code	CSB-PA070249
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TD19
Immunogen	Synthesized peptide derived from Human Nek9 around the non-phosphorylation site of T210.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NEK9; KIAA1995; NEK8; NERCC; Serine/threonine-protein kinase Nek9; Nercc1 kinase; Never in mitosis A-related kinase 9; NimA-related protein kinase 9; NimA-related kinase 8; Nek8
Immunogen Species	Homo sapiens (Human)
Target Names	NEK9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.