



GCNT3 Antibody

Product Code	CSB-PA070267
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95395
Immunogen	Synthesized peptide derived from the Internal region of Human GCNT3.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GCNT3; Beta-1; 3-galactosyl-O-glycosyl-glycoprotein beta-1,6-N-acetylglucosaminyltransferase 3; C2GnT-mucin type; C2GnT-M; hC2GnT-M; Core 2/core 4 beta-1,6-N-acetylglucosaminyltransferase; C2/4GnT
Immunogen Species	Homo sapiens (Human)
Target Names	GCNT3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.