




MAPKAPK2 Antibody

Product Code	CSB-PA070278
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49137
Immunogen	Synthesized peptide derived from Human MAPKAPK-2 around the non-phosphorylation site of T334.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MAPKAPK2; MAP kinase-activated protein kinase 2; MAPK-activated protein kinase 2; MAPKAP kinase 2; MAPKAP-K2; MAPKAPK-2; MK-2; MK2
Immunogen Species	Homo sapiens (Human)
Target Names	MAPKAPK2

Image



Western blot analysis of 293T 3T3 lysis using MAPKAPK-2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.